User's Manual Blood Glucose Meter

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© Contents

Introduction	3
Chapter 1: Your New System	5
Chapter 2: Control Tests	13
Chapter 3: Blood Glucose Tests	19
Chapter 4: Meter Settings	29
Chapter 5: Review Your Data	41
Chapter 6: Wireless Communication and Meter Pairing	47
Chapter 7: Data Transfer and PC Reports	61
Chapter 8: Meter Maintenance and Troubleshooting	77
Chapter 9: Technical Information	87
Chapter 10: Guarantee	93
Index	95

The Accu-Chek Aviva Connect System

Intended use

The Accu-Chek Aviva Connect meter is designed to be used with the Accu-Chek Aviva test strip to quantitatively measure glucose in fresh capillary whole blood as an aid in monitoring the effectiveness of glucose control.

Capillary whole blood for testing of blood glucose can be obtained from fingertip and approved alternate sites (for example, forearm). Refer to the Alternate Site Testing (AST) section of this manual for approved alternate sites and associated limitations. The Accu-Chek Aviva Connect meter with the Accu-Chek Aviva test strips provide a complete test system that is meant for in vitro diagnostic use by a single person and should not be shared. The system is not for use in diagnosis or screening of diabetes mellitus or for neonatal use.

For use only with the Accu-Chek Aviva test strips and control solutions

Suitable for self-testing

The system includes:

Accu-Chek Aviva Connect meter with batteries, Accu-Chek Aviva test strips*, Accu-Chek Aviva control solutions*, and Accu-Chek Connect diabetes management app.**

*Some items may not be included in the kit. They are a separate purchase.

** Download and install the app from the app store.

A WARNING

- Choking hazard. Small parts. Keep away from children under the age of 3 years.
- Any object coming into contact with human blood is a potential source of infection (see: Clinical and Laboratory Standards Institute: Protection of Laboratory Workers from Occupationally Acquired Infections; Approved Guideline – Third Edition; CLSI document M29-A3, 2005).

🕂 WARNING

The lancing device is for single patient use only and should not be shared.

Why Regular Blood Glucose Testing Is Important

Testing your blood glucose regularly can make a big difference in how you manage your diabetes every day. We have made it as simple as possible.

Important Information About Your New Meter

- The meter prompts you to choose a language the first time you turn it on.
- Check the time and date on the meter before you begin testing. Adjust the time and date if necessary.
- Sample data displays are shown throughout the manual. Your data will differ.
- If you follow the steps in this manual but still have symptoms that do not seem to match your test results, or if you have questions, talk to your healthcare professional.

The Accu-Chek Aviva Connect Meter

Front View

1. Back Button

Returns to a previous display or field.

 Up Arrow and Down Arrow Button Press to move between menu options or to increase or decrease numbers.

3. Display

Shows results, messages, and test results stored in memory.

4. Power/Set/OK Button

Turns meter on or off and sets options.

5. Test Strip Slot

Insert test strip here.

Back View

6. Battery Drawer Pull out to replace batteries.

Side View

7. Micro USB Port Transfers data from the meter to a computer (PC).

Your New System

10. Yellow Window Touch blood drop or control solution here. Connects the meter to a PC.

*Some items may not be included in the kit. They are a separate purchase.

1

Button Functions

Here are the functions of the back, arrow, and Power/Set/OK buttons on the meter. These functions are used throughout this manual. See the chapter **Meter Settings** for specific instructions on setting up the meter.

Button	Function
Back Button	Return to the previous display. Return to the previous field.
Up Arrow and Down Arrow Button	Navigate up and down in a menu. Increase or decrease a number.
Power/Set/OK Button	Press briefly to turn the meter on. Press and hold to turn the meter off. Press to select an option. Press to move to the next field or display. Press to save an option. With the meter off, press and hold to check the meter display.
	Display Check Display Check

Meter Menus

1

- Title of display or menu
- Menu options
- Scroll down
- Scroll up

Display	Description
	There are more menu options listed beneath the last option. Press \blacksquare on the meter to view the options.
	There are more menu options listed above the first option. Press \blacksquare on the meter to view the options.
	There are more menu options listed above and below the options. Press \blacksquare or \blacksquare on the meter to view the options.

Your New System 1

	1	i
E	1	l
	I	I
	I	L
		L

Display	Description
9:38am 12/11/12 Main Menu Test My Data	Highlighted option (Test) Press 🚾 to enter the Test menu.
Time/Date HH MM A/P 9:38 am MM DD YY 12/11/12	Highlighted field (Hour) Press or or to increase or decrease the hour. Press

Icons

1

Here is a list of the icons on the display.

lcon	Description
1	Above target range
Ì	After meal
⊨ ≕ı	Bedtime
Ù	Before meal
Ţ	Below target range
•	Blood glucose test
\checkmark	Checkmark
C	Control bottle
×	Control test not OK
\checkmark	Control test OK
(\mathbf{X})	Error

Your New System 1

lcon	Description
\otimes	Fasting
	Flight mode
?	Help
Ē	Low battery
1	My data
	No comment
*	Other
\odot	Overall
~	Selected setting
	Settings
	Temperature warning
\triangle	Warning
\$	Within target range

Set the Language

The meter prompts you to choose a language the first time you turn it on.

NOTE

If you select the wrong language and cannot correct it, contact Accu-Chek Customer Care at 1-800-363-7949.

When to Perform a Control Test

Performing a control test lets you know the meter and test strips are working properly. You should perform a control test when:

- you open a new test strip box.
- you left the test strip container open.
- you think the test strips are damaged.
- you want to check the meter and test strips.
- the test strips were stored in extreme temperatures, humidity, or both.
- you dropped the meter.
- your test result does not match how you feel.
- · you want to check if you are performing the test correctly.

About the Control Solutions

- Use only Accu-Chek Aviva control solutions.
- Close the control solution bottle tightly after use.
- Write the date you open the control solution bottle on the bottle label. The control solution must be discarded 3 months from the date the control solution bottle was opened (discard date) or on the use by date on the bottle label, whichever comes first.
- Do not use control solution that is past the use by or discard date.
- Refer to the control solution package insert for control solution storage conditions.
- The meter automatically recognizes the difference between the control solution and blood.
- The control results are not displayed in memory.
- The control solution can stain fabric. Remove stains by washing with soap and water.

Performing a Control Test

You need the meter, a test strip, and control solution Level 1 or Level 2.

Check the use by date on the test strip container. Do not use test strips past the use by date.

2

Insert the test strip into the meter in the direction of the arrows. Place the meter on a flat surface.

The meter turns on. **Preparing to test** appears.

Apply drop appears.

2

Select the control solution to test. You will enter the level later in the test.

Remove the bottle cap. Wipe the tip of the bottle with a tissue. Squeeze the bottle until a tiny drop forms at the tip.

5

Touch the drop to the **front edge** of the yellow window of the test strip. Do not put control solution on top of the test strip.

Analyzing appears when there is enough control solution in the test strip.

2 Control Tests

Control Result and the control bottle icon appear.

Press \square or \square to select the control level you tested. If you do not select a level, the control level is saved as --.

If any other action is taken at this time, the control result is saved with no control level.

Press OK.

✓ and Within range appear if the control result is within range.

X and **Out of range** appear if the control result is out of range.

2

8	9
Wipe the tip of the bottle with a tissue. Cap the bottle tightly.	Remove and discard the used test strip.

NOTE

The meter turns off 90 seconds after a successful test or 15 seconds after the test strip is removed, provided no other action is taken.

Understanding Out-of-Range Control Results

A WARNING

The control ranges are printed on the test strip container label. If the control result is out of range, check this list to help solve the problem.

Troubleshooting Checks	Action
Were the test strips or control solutions expired?	Discard the test strips or control solution if either is past the use by date. If the control solution was opened more than 3 months ago, discard it. Repeat the control test with an unexpired test strip and an unexpired control solution.
Did you wipe the tip of the control solution bottle before use?	Wipe the tip of the bottle with a tissue. Repeat the control test with a new test strip and a fresh drop of control solution.
Were the caps on the test strip container and the control solution bottle always closed tightly?	Replace the test strips or control solution if you think either was uncapped for some time. Repeat the control test.
Was the test strip used immediately after it was removed from the test strip container?	Repeat the control test with a new test strip and a fresh drop of control solution.
Were the test strips and control solutions stored in a cool, dry place?	Repeat the control test with a properly stored test strip or control solution.
Did you follow the directions?	Read the chapter Control Tests and repeat the control test.
Did you choose the correct control solution level, either 1 or 2, when you performed the control test?	If you chose the wrong control solution level, you can still compare the control result to the range printed on the test strip container.
Are you still unsure of the problem?	Contact Accu-Chek Customer Care at 1-800-363-7949.

3

\land WARNING

Blood glucose results can be displayed in either mg/dL or mmol/L. The back label of the meter shows the unit of measurement. Meters for use in Canada have mmol/L as the unit of measurement. If the meter shows the wrong unit, contact Accu-Chek Customer Care at 1-800-363-7949. If you do not know which unit of measurement is correct for you, contact your healthcare professional. Using the wrong unit of measurement may cause misinterpretation of your actual blood glucose level and may lead to improper therapy.

Using the Accu-Chek Aviva Connect System

- Use only Accu-Chek Aviva test strips.
- Use the test strip immediately after removing it from the test strip container.
- . Do not apply blood or control solution to the test strip before inserting it into the meter.
- Close the test strip container tightly immediately after removing a test strip to protect the test strips from humidity.
- Store the unused test strips in their original test strip container with the cap closed.
- Check the use by date on the test strip container. Do not use the test strips after that date.
- Store the test strip container and meter in a cool, dry place such as a bedroom.
- · Refer to the test strip package insert for test strip storage and system operating conditions.

\land WARNING

Do not store test strips in high heat and moisture areas (bathroom or kitchen)! Heat and moisture can damage test strips.

Performing a Blood Glucose Test with Blood from Your Fingertip

NOTE

- A blood glucose test cannot be performed while the meter is connected to a PC with a USB cable.
- Before you perform your first blood glucose test, set up the meter correctly.
- You need the meter, a test strip, and a lancing device with a lancet inserted to perform a blood test.
- There are 2 ways to start a blood glucose test.
 - Insert a test strip into the meter.
 - Turn the meter on by briefly pressing \mathbf{M} . Select **Test** > \mathbf{M} .

Wash your hands with warm soapy water and dry thoroughly.

Prepare the lancing device.

Check the use by date on the test strip container.

Do not use test strips past the use by date.

Insert the test strip into the meter in the direction of the arrows.

The meter turns on. **Preparing** to test appears.

When **Apply drop** appears, perform a fingerstick with the lancing device.

Gently squeeze your finger to assist the blood flow. This helps you get a blood drop.

Meter Settings for details on

Touch the end of the test strip to the blood drop. Do not put blood on top of the test strip.

Analyzing appears when there is enough blood in the test strip.

Last Result.

Add Comment appears. Press or to highlight the comment. Press ox.

The final result appears.

With the comment highlighted, either press I™ to set the comment or press I™ to return to Add Comment to change the comment. If the result field is highlighted, press I™ to set the option and return to Main Menu. 3

Here is a list of icons that can be added to a blood glucose result.

lcon	Description	
Ŭ	Before meal	
Ť	After meal	
\otimes	Fasting	
	Bedtime	
*	Other	
	No entry	1. You do not want to add a comment.
		2. You want to remove a comment for the current blood glucose result.

9

Remove and discard the used test strip.

If you see either of these results:

9:38am 12/11/12 LO Add Comment 7:38am 4.7mmol/L	Blood glucose may be lower than the measurement range of the system. See Unusual Blood Glucose Results in this chapter.
9:38am 12/11/12 HII Add Comment 7:38am 4.7mmol/L	Blood glucose may be higher than the measurement range of the system. See Unusual Blood Glucose Results in this chapter.

Performing a Blood Glucose Test with Blood from Your Palm, Forearm, or Upper Arm (Alternate Site Testing)

You have the option of obtaining a blood sample from other sites on your body besides the fingertip. Alternate sites include the palm, forearm, and upper arm.

Blood obtained from the fingertip can be used at any time to perform a blood glucose test.

If blood from the palm, forearm, or upper arm is used, there are certain times when testing is not appropriate. This is because your blood glucose level changes faster in your fingertip than in the palm, forearm, and upper arm. These differences may cause you to misinterpret your actual blood glucose level, leading to improper therapy and potential adverse health effects.

Read the next section before you try testing from the palm, forearm, or upper arm.

You may perform a palm, forearm, or upper arm test	 immediately before a meal. while fasting.
You may NOT perform a palm, forearm, or upper arm test	 up to 2 hours following a meal, when blood glucose values can rise quickly. after injecting bolus insulin, when blood glucose values can decrease rapidly. after exercise. if you are sick. if you are sick. if you think your blood glucose is low (hypoglycemia). if you sometimes do not notice when your blood glucose is low.

If you are interested in AST, talk to your healthcare professional first.

To obtain an AST cap and detailed AST instructions, contact Accu-Chek Customer Care at 1-800-363-7949.

Unusual Blood Glucose Results

If your blood glucose result does not match how you feel, check this list to help solve the problem.

Troubleshooting Checks	Action
Were the test strips expired?	Discard the test strips if they are past the use by date. Repeat the blood glucose test with an unexpired test strip.
Was the cap on the test strip container always closed tightly?	Replace the test strips if you think the test strip container was uncapped for some time. Repeat the blood glucose test.
Was the test strip used immediately after it was removed from the test strip container?	Repeat the blood glucose test with a new test strip.
Were the test strips stored in a cool, dry place?	Repeat the blood glucose test with a properly stored test strip.
Did you follow the directions?	See the chapter Blood Glucose Tests and repeat the blood glucose test. Contact Accu-Chek Customer Care at 1-800-363-7949 if you still have problems.
Are the meter and test strips working properly?	Perform a control test. See the chapter Control Tests for instructions.
Are you still unsure of the problem?	Contact Accu-Chek Customer Care at 1-800-363-7949.

Symptoms of Low or High Blood Glucose

Being aware of the symptoms of low or high blood glucose can help you understand your test results and decide what to do if they seem unusual.

Low blood glucose (hypoglycemia): Symptoms of hypoglycemia may include, but are not limited to, anxiety, shakiness, sweating, headache, increased hunger, dizziness, pale skin color, sudden change in mood or irritability, fatigue, difficulty concentrating, clumsiness, palpitations, and/or confusion.

High blood glucose (hyperglycemia): Symptoms of hyperglycemia may include, but are not limited to, increased thirst, frequent urination, blurred vision, drowsiness, and/or unexplained weight loss.

\land WARNING

If you are experiencing any of these symptoms, or other unusual symptoms, test your blood glucose from the fingertip. If your blood glucose result is displayed as LO or HI, contact your healthcare professional immediately.

4

Overview

You can adjust the following settings in the meter for your personal preferences. Refer to the sections later in this chapter for details and how to set the options.

Setting	Options	Function
Time/Date	Time / Date	Set the time and date.
Beeper	On / Off	Select On or Off.
USB Default	Data Transfer / PC Reports	Select the desired activity for the data in the meter when the USB cable is plugged into the meter and a PC.
		Data Transfer – copies the data in the meter to a compatible software program.
		PC Reports – the meter generates data reports that open in an Internet browser and can be printed.
Target Range	On / Off	Set Markers
		Select On or Off.
		$\mbox{On}-\mbox{blood}$ glucose results are marked as above, within, or below based on the target range set in the meter.
		\ensuremath{Off} – no target range markers appear with a blood glucose result.
	3,9–8,9 mmol/L (pre-set target range)	Set Range
		Select the blood glucose target range appropriate for you. Consult your healthcare professional for the appropriate target range for you.

Meter Settings

Setting	Options	Function
Last Result	On / Off	Select whether the previous blood glucose result (within the past 24 hours) ③ appears with the current blood glucose result ①. On – the previous blood glucose result appears with the current blood glucose result.
		Off - only the current blood glucose result appears.
		9:38am 12/11/12 5.8 mmol/L Add Comment 7:38am 4.7mmol/L 2
Language		Select the language for the meter.
Wireless		Select wireless communication settings. See the chapter Wireless Communication and Meter Pairing.

Δ

Time/Date (Main Menu > Settings > Time/Date)

4

Press of to save and return to **Settings**.

Beeper (Main Menu > Settings > Beeper)

The beeper prompts you:

4

- when a test strip is inserted.
- to apply blood or control solution to the test strip.
- when enough blood or control solution is drawn into the test strip.
- when the blood glucose or control test is complete.
- when the meter is turned on.
- when a button is pressed.
- when the batteries are inserted.
- when there are no stored blood glucose results or errors or there is an invalid record.
- if an error occurred (even if the beeper is off, it still beeps for an error).

Turn the meter on by briefly pressing . From Main Menu, press a to highlight Settings. Press .

Press or to highlight **On** or **Off**. Press to move \checkmark to the option. Press **C** to set the option and return to **Settings**. Δ

USB Default (Main Menu > Settings > USB Default)

The USB Default option determines how the data in the meter is handled whenever you plug the USB cable into the meter and a PC. Your default option remains the same until you change it.

There are 2 options to display and analyze results on a PC.

- 1. Data Transfer the data is read by special software for diabetes management on the computer.
- 2. PC Reports the meter generates data reports that open in an Internet browser and can be printed.

See the chapter Data Transfer and PC Reports for detailed information on data management.

Δ

move ✓ to the option. Press III to set the option and return to Settings.

Target Range (Main Menu > Settings > Target Range > Set Markers / Set Range)

Set Markers

- If Set Markers is off, the markers do not appear with the blood glucose results.
- Any marker remains with the blood glucose result in the Logbook, even if you later turn Set Markers off.
- If Set Markers is on, the following markers appear with blood glucose results.

Marker	Meaning
Ţ	The blood glucose result is below the target range.
\$	The blood glucose result is within the target range.
1	The blood glucose result is above the target range.

Meter Settings

Set Range

- Consult your healthcare professional for the appropriate target range for you.
- The target range can be set from a lower limit of 2,8-5,5 mmol/L to an upper limit of 5,6-11,1 mmol/L.

\land WARNING

Δ

This function is no substitute for hypoglycemia training by your healthcare professional.

Press or to highlight On or Off. Press to move \checkmark to the option. Press to set the option and return to Target Range.

Press To highlight **Set Range**. Press **M**.

Press or to adjust the lower limit of the range. Press

Press or to adjust the upper limit of the range. Press to set the option. Target Range appears. Press to return to Settings.

Δ

Last Result (Main Menu > Settings > Last Result)

Select whether the previous blood glucose result appears with the current blood glucose result. Test results older than 24 hours do not appear.

- On the previous blood glucose result appears with the current blood glucose result.
- Off only the current blood glucose result appears.

Δ

Language (Main Menu > Settings > Language)

Choose the language that appears on the meter.

Overview

- Blood glucose results are stored from the newest to the oldest.
- The meter automatically stores up to 750 blood glucose results in memory with the time and date of the test and any test result markers.
- Once more than 750 blood glucose results are in memory, adding a new blood glucose result deletes the oldest blood glucose result.
- Only test results that have been assigned a fasting, before meal, after meal, or bedtime marker are included in the average for that marker.
- All test results are included in the overall 7-, 14-, 30-, and 90-day averages regardless of what comment is assigned.
- Control results are stored in memory but cannot be reviewed on the meter. To view stored control results, transfer them to a compatible software application.
- Control results are not included in the averages or blood glucose reports.
- Once more than 30 control results are in memory, adding a new control result deletes the oldest control
 result.

\land WARNING

Do not change your therapy based on an individual test result in memory. Talk to your healthcare professional before changing therapy based on test results in memory.

Logbook (Main Menu > My Data > Logbook)

Press or to scroll through **Logbook**.

most recent result

2nd most recent result

3rd most recent result

4th most recent result

Averages (Main Menu > My Data > Averages)

Transfer Data (Main Menu > My Data > Data Transfer/PC Reports)

This menu allows you to choose either **Data Transfer** or **PC Reports**, regardless of what the **USB Default** setting is or if you do not set the **USB Default**. For example, if the **USB Default** is **PC Reports**, but your healthcare professional wants to transfer your data, navigate to this menu and select **Data Transfer**. If the **USB Default** is not set, you can use this menu to select where to view your data each time.

3

Connect the USB cable to the meter.

Overview

You can wirelessly and automatically synchronize your diabetes information with a mobile device and your personal ACCU-CHEK Connect online account. The process of creating a connection between the meter and the other device is called pairing. You need an application that can accept the meter's data.

Setting	Options	Function
Flight Mode	On / Off	Select whether wireless communication is available.
		On – wireless communication is not available.
		Off - wireless communication is available.
Default Device	List of devices paired with the meter	If more than 1 device is paired, select the device with which the Auto-Send and Sync Time features will communicate.
Auto-Send	On / Off	Select whether data is automatically sent to the default paired device after each test.
		$\mbox{On}-\mbox{data}$ is automatically sent to the default paired device.
		Off – data is not automatically sent to the default paired device.
Sync Time	On / Off	Select whether to synchronize the time and date to the default paired device.
		\mbox{On} – the time and date on the meter synchronize to the time and date on the default paired device.
		Off – the time and date on the meter do not synchronize to the time and date on the default paired device.
Pairing	Pair Device / Delete Pairing	Select whether to pair a device or to delete a paired device.

Wireless (Main Menu > Settings > Wireless)

The steps below are for first time pairing.

6

Up to 5 devices can be paired with the meter at any time.

Each device must be paired with the meter one at a time.

The meter and the device to be paired should be within 1 meter of each other.

Wireless Communication and Meter Pairing

The other device displays a list of found devices, including the meter. Select Accu-Chek from the list. When prompted, enter the meter code from Step 5 into the device to be paired.

6

When pairing is complete, a checkmark and the name of the paired device appear. Press at to return to **Wireless**.

Flight Mode (Main Menu > Settings > Wireless > Flight Mode)

Select whether wireless communication is available or not. When **Flight Mode** is on, the airplane icon appears in the title bar and wireless communication is not available.

Press or to highlight **On** or **Off**. Press to move to the option. Press to set the option and return to **Wireless**.

If more than 1 device is paired, select the default paired device for the Auto-Send and Sync Time features.

device you want as the default device. Press or to set the option and return to **Wireless**.

6

Auto-Send (Main Menu > Settings > Wireless > Auto-Send)

Select whether data is automatically sent to the default paired device after each test.

Press or to highlight **On** or **Off**. Press to move \checkmark to the option. Press **C** to set the option and return to **Wireless**.

Sync Time (Main Menu > Settings > Wireless > Sync Time)

Select whether to synchronize the time and date to the default paired device.

the option. Press **to** set the option and return to **Wireless**.

6

Pairing Additional Devices (Main Menu > Settings > Wireless > Pairing > Pair Device)

This procedure is to add another paired device (assumes there is already 1 paired device).

You can pair 5 devices. If you have 5 paired devices, you must delete a pairing before you can add another device.

The device you choose remains as the default until you change it.

6

Wireless Communication and Meter Pairing

Pairing
Pair Device
Delete Pairing
⑦ Help

Press or to highlight **Pair Device**. Press **C**.

If there are less than 5 paired devices, **Pair Device** appears. Go to Step 5.

If there are already 5 paired devices, **Pair Device** appears. You must delete a paired device before you can proceed to Step 5. Select **Yes**. Press or to highlight the device to delete. Press **C**. Proceed to Step 5.

5

4

Prepare the other device for pairing. Refer to the other device's instructions for information about pairing.

The meter displays its code.

7

The other device displays a list of found devices, including the meter. Select Accu-Chek from the list. When prompted, enter the meter code from Step 6 into the device to be paired.

Wireless Communication and Meter Pairing

6

When pairing is complete, \checkmark and the name of the paired device appear. Press $\square K$.

If there is more than 1 paired device, you will be asked if you want the new device to be the default device. Press or to highlight No or Yes. Press or to set the option and return to Pairing. Press to return to Wireless. This procedure is to delete a paired device (assumes there is at least 1 paired device).

You can pair 5 devices. If you have 5 paired devices, you must delete a pairing before you pair to another device.

The device you choose remains as the default until you change it.

6

6 Wireless Communication and Meter Pairing

4	5
Pairing	Delete Pairing
Pair Device	NOKIA1234
Delete Pairing	NOKIA2345
⑦ Help	NOKIA3456
	NOKIA4567
Press ➡ or ➡ to highlight Delete Pairing. Press .	Press or to highlight the device to delete. Press of to set the option

Delete Pairing

Press OK to

choose new default device.

6

If the selected device is not the default device, the pairing is deleted. \checkmark and the name of the deleted device appear.

If the selected device is the default device, you are asked to select a new default device. If there is only 1 device remaining, it becomes the default device. The list of paired devices appears.

6

Highlight the new default device and press \square . \checkmark and the name of the deleted device appear, along with the new default device.

6 Wireless Communication and Meter Pairing

You have 2 options to display and analyze blood glucose results on a PC.

Data Transfer - this option transfers the data to special software for diabetes management in a PC.

PC Reports - the meter generates data reports that open in an Internet browser and can be printed.

The option you select remains as the default until you change it.

NOTE

- You cannot perform a blood glucose test while the meter is connected to a computer with a USB cable. If you connect the meter to a PC when a test is in progress, the test is cancelled.
- The meter has a port on the side for the small end of the USB cable. The large end of the USB cable is inserted into the USB port of a PC.

Set Data Transfer as the Default (Main Menu > Settings > USB Default > Data Transfer)

Set PC Reports as the Default (Main Menu > Settings > USB Default > PC Reports)

View the Data in Compatible Diabetes Management Software

These steps assume that the USB default option is set to Data Transfer.

The meter can be off or on.

The meter transfers the data to the software.

NOTE

If you would like to view the data in PC Reports instead, follow these steps:

- 1. Unplug the USB cable. Connection Lost appears.
- 2. Press 🖪 to return to Main Menu.
- 3. Select My Data>PC Reports.
- 4. Reconnect the USB cable.
- 5. PC Reports appears.

View the Data on a PC (in an Internet Browser)

These steps assume that the USB default option is set to PC Reports.

The meter can be off or on.

Once the meter is disconnected from the PC, the data disappears from the PC (but not the meter) unless you save it on the PC.

Plug the small end of the USB cable into the meter.

Plug the large end of the USB cable into a USB port on the PC.

If the meter is off, it turns on and displays **PC Reports**.

3	4	5
Open the file manager (for example, Microsoft Windows Explorer) on the PC. The meter appears as a drive (USB storage device) in the file manager.	Double click on the ACCU-CHEK drive icon.	Click on Start.html . The Internet browser opens and the default reports appear.

NOTE

If you would like to view the data in special diabetes management software instead, follow these steps:

- 1. Unplug the USB cable. Connection Lost appears.
- 2. Press 🖾 to return to Main Menu.
- 3. Select My Data>Data Transfer>USB Cable.
- 4. Reconnect the USB cable.
- 5. Data Transfer appears.

Shortcut

Here is a shortcut to get to **Data Transfer** or **PC Reports**. This is a quick way to transfer the data to software or to view the meter's data on a PC.

- 1. Turn the meter off.
- 2. Press and hold both 🗳 and 🖾 until Activity appears.
- 3. Select either Data Transfer or PC Reports.

View PC Reports

Overview

- An Internet connection is not necessary to view PC Reports.
- When the meter and PC successfully connect, the following reports are displayed in the Internet browser:
 - Trend Report This report shows the trend for several test results over the selected time period (past 3, 7, 14, 30, or 90 days).
 - Standard Day Report This report shows all data in a 24-hour grid.
 - Standard Week Report This report shows all blood glucose results according to the time when the test
 was performed and the day of the week.
 - List Report The list report (record list) shows the test results sorted by date and time of the test.

Data Transfer and PC Reports

Trend Report

27.04.2013-10.05.2013

Statistics

Below the chart of a report, you will find a statistical analysis of all test results plotted with the following information:

- Number of tests
- Average
- Standard deviation The standard deviation is the variance of the analyzed results.
- Highest value
- · Lowest value
- Average tests per day Average number of blood glucose tests per day
- HI or LO Blood glucose results outside the measuring range
- Low blood glucose index or High blood glucose index Further information can be found in References
- Above target range Blood glucose results above the target range
- Within target range Blood glucose results within the target range
- Below target range Blood glucose results below the target range

Key

х	Blood glucose result
Х	Several blood glucose results
	Average blood glucose results in the selected time period
^	Blood glucose result above the chart range
^ (red caret)	Blood glucose may be higher than the measurement range of the system
V	Blood glucose may be lower than the measurement range of the system
(green)	Your personal blood glucose target range (shown as a green bar on the PC monitor)
(black)	Non-work days

Print Report

Do not use the print function of the Internet browser. Use the **Print reports** button instead.

Excluded Data

The following blood glucose results are not included in a report:

- · Test results outside the selected time period
- · Control results
- · L0 or HI test results
Trend Report

The Trend Report shows the trend of blood glucose results over the selected time period.

The date appears on the horizontal x-axis. The blood glucose results appear on the vertical y-axis. The test results are connected by a thin black line in chronological order.

If you set a target range in the meter, it appears as a green bar on the chart. Non-work days are marked with diagonal slashes on the horizontal x-axis. The trend of the day-to-day average blood glucose result is represented by a thick black line.

Standard Day Report

The Standard Day Report makes it easier to recognize daily patterns. All blood glucose results are placed on a 24-hour grid. Therefore, all tests performed at (approximately) the same time are shown at the same position on the horizontal time axis.

Blood glucose results are connected by a thin black line in chronological order. A thick black line represents the trend of the average level (in intervals of 1 hour if a test result falls in each interval).

Standard Week Report

The Standard Week Report makes it easier to recognize weekly patterns. For example, you might find trends brought about by your occupation.

All blood glucose results are plotted on the chart according to the time and day the test was performed. Blood glucose results are connected by a thin black line in chronological order. A thick black line represents the trend of the average level for each day.

The number of tests and the daily blood glucose average are listed below the chart.

List Report

List Report

08.05.2013-10.05.2013

Date and Time		Blood Glucose (mmol/L)	Events
05/10/2013	07:30am	11.7	Before meal
	01:01am	9.9	
05/09/2013	08:00pm	9.4	After meal
	05:45pm	3.9	Before meal
	03:00pm	3.6	After meal
	12:27pm	2.2	Before meal
	10:18am	16.6	After meal
	07:37am	17.2	Before meal
	01:10am	10.4	
05/08/2013	10:00pm	3.9	After meal
	08:10pm	3.9	After meal
	05:27pm	4,7	Before meal
	03:17pm	3.7	After meal
	12:25pm	1.8	Before meal
	10:01am	16.6	After meal

The List Report shows the test results sorted by date and time of the test. All blood glucose results are listed chronologically with any additional information about the test result.

The list contains the following columns:

- Date and Time
- Blood Glucose
- Events Event connected to this test result

Analyzing Data in External Applications

If you want to analyze the test results using external software, you can save the data as a CSV file (Comma Separated Values). CSV files can be opened with a text editor or spreadsheet program.

The CSV file contains all test results saved in the meter. Test results transferred at an earlier date are transferred again.

1. Press the **Save file** button ӣ on the user interface.

Depending on the configuration of the PC, the CSV file may be directly opened in a spreadsheet program. In this case, you can save the data using the **Save** function of the spreadsheet program.

 In the dialog box that opens, select the option to save the file. Where you find the CSV file on the PC depends on the settings in the operating system for data download.

The CSV file contains the following information:

- Serial number: Serial number of the meter.
- Download date, download time: date and time when the meter transferred the test results to the computer.
- Date, time, result and unit of the results saved in the meter.
- Flags added to the results, indicated by an X.

In the CSV file, the date is always displayed as DD.MM.YYYY and the time as 24-hour format (hh:mm). The time format set in the meter has no influence on the format of the date and time in the CSV file.

Security Settings within the Internet Browser

The settings of the Internet browser can influence working with reports. The reports use pages with active content (JavaScript). This active content can be suppressed by security settings in the browser, causing warnings or restricted functionality. If this happens, check the Internet browser settings.

In many cases, you can create different security settings for using the Internet and working with reports (for example, at user login to the PC or by defining user profiles in the browser).

If you select the Internet browser security settings appropriately (for example, **Allow active content to run in files on My Computer**), you can work with reports without any restrictions.

NOTE

Some newer Internet browsers do not allow the **Save** function (e.g., Mozilla Firefox version 15 and higher). The files can be saved from Windows Explorer and subsequently opened using programs that are compatible with *.CSV files.

Troubleshooting

Troubleshooting Check	Action
The ACCU-CHEK drive symbol with the start.html file does not appear on the PC.	Check whether PC Reports is selected as the default.
	Check whether the PC or operating system supports data transfer via USB.
	Check whether the USB connector is firmly plugged into the correct port on the PC.
The meter is still not detected as a drive.	Plug the meter into a different USB port on the PC.

Meter Maintenance

The meter automatically tests its own systems every time you turn it on and lets you know if something is wrong. See **Error Messages** in this chapter.

If you drop the meter or think the results are not accurate, contact Accu-Chek Customer Care at 1-800-363-7949.

Changing the Batteries

Use your thumb to slide the battery drawer out of the meter.

Remove the old batteries and place the new ones in the battery drawer with the (+) side facing down.

8

Slide the battery drawer back into position until it locks into place.

NOTE

- The meter uses two 3-volt lithium batteries, coin cell type CR2032. This type of battery can be found in many stores. It is a good idea to have spare batteries available.
- Always replace both batteries at the same time and with the same brand.
- The logbook data is saved when you replace the batteries.

Meter Cleaning and Disinfecting

8

What is the difference between cleaning and disinfecting?

Cleaning is the removal of dirt from the meter.

Disinfecting is the removal of most, but not all, disease-causing and other types of microorganisms (bloodborne pathogens) from the meter.

Approved Cleaning and Disinfecting Product

The following active ingredients have been approved for cleaning and disinfecting the meter:

 A mixture of alkyl dimethyl ethylbenzyl ammonium chloride (0.125 %), benzalkonium chloride (0.125 %), and isopropyl alcohol (14.85 %).

We recommend Sani-Cloth Plus Germicidal Disposable Cloths. This product contains the approved active ingredients and is licensed for use in Canada.

Roche has tested the approved product for a total of 520 cleaning and disinfecting cycles (260 cleaning and 260 disinfection cycles), which is equal to cleaning and disinfecting once per week for a 5-year period.

For more information see:

- FDA Public Health Notification: Use of Fingerstick Devices on More than One Person Poses Risk for Transmitting Bloodborne Pathogens; Initial Communication, (2010). http://www.fda.gov/MedicalDevices/Safety/AlertsandNotices/ucm224025.htm Accessed March 26, 2014.
- 2. CDC website on "Infection Prevention during Blood Glucose Monitoring and Insulin Administration". http://www.cdc.gov/injectionsafety/blood-glucose-monitoring.html Accessed March 26, 2014.

NOTE

- Always wash hands thoroughly with soap and water before handling the meter, lancing device, or test strips.
- For technical assistance or questions on cleaning and disinfecting, contact Accu-Chek Customer Care at 1-800-363-7949.

Cleaning and Disinfecting the Meter

To clean and disinfect without damaging the meter, follow these procedures carefully.

When to Clean and Disinfect the Meter

- Clean the meter to remove visible dirt or other material prior to disinfecting.
- Clean and disinfect the meter at least once per week and when blood is present on the surface of the meter.
- Clean and disinfect the meter before allowing anyone else to handle the meter. Do not allow anyone else to use the meter on themselves for testing purposes.

What to Clean and Disinfect

The following parts of the meter should be cleaned and disinfected:

- The area around slots and openings (do not get any moisture in slots or openings)
- The meter display
- The entire meter surface

How to Clean and Disinfect the Meter

A WARNING

Failure to follow these instructions will damage the meter and stop it from working properly.

- DO NOT clean or disinfect the meter while performing a blood glucose or control test.
- D0 NOT get any moisture in slots or openings.
- DO NOT spray anything onto the meter.
- DO NOT immerse the meter in liquid.
- Always use the same product for both cleaning and disinfecting.

Wash hands thoroughly with soap and water.

Turn the meter off and wipe the entire meter surface with a Sani-Cloth Plus. Carefully wipe around the test strip slot and other openings.

Make sure that no liquid enters any slot or opening.

3

A separate Sani-Cloth Plus should be used for cleaning and disinfection. For disinfecting the meter, get a new cloth and repeat step 2, making sure that the surface stays wet for 3 minutes.

Make sure that no liquid enters any slot or opening. Ensure the meter is completely dry before performing a blood glucose test.

8

Wash hands thoroughly with soap and water.

Error Messages

A WARNING

- Never make therapy decisions based on an error message.
- If you have any concerns or see any other error message, contact Accu-Chek Customer Care at 1-800-363-7949.

The meter will not turn on or the display is blank.	▲Test Not Allowed	▲Connection Lost
Batteries are dead.	Remove strip	Reconnect USB
Insert new batteries.	or LISB cable	cable and retry.
Display is damaged. / Meter is defective.		
Contact Accu-Chek Customer Care at 1-800-363-7949.	The meter is connected to a PC and a test cannot be performed.	The connection between the meter and PC was lost.
Extreme temperatures.	EITHER remove the USB cable	Disconnect and reconnect the
Move the meter to a more temperate area.	the test strip and start a data transfer.	connection. Contact Accu-Chek Customer Care at 1-800-363-7949 if the connection is lost again.

A No Connection Connect USB cable to continue.	▲Transfer Failed Ensure paired device is within range and turned on. Blood glucose results were not transferred to a paired device.	Connection Failed Check PC or USB cable. Data could not be transferred from the meter to the PC.
Connect the USB cable and retry the connection.	Make sure the paired device is within range of the meter and turned on.	Check the PC or USB cable.

A Transfer Not Allowed Wireless communication not allowed in flight mode.	A Results Excluded One or more results may be excluded from averages.	▲ Invalid Date Check the date.
Data cannot be sent to a paired	One or more blood glucose results are excluded from the selected averages because the	The date entered is not valid.
device because the meter is in Flight Mode.		Enter the correct date.
Retry the data transfer when the meter is not in Flight Mode .	results are invalid or out of the system measurement range.	

A Flight Mode Active Wireless is off. Selected setting was saved.	A Flight Mode Active Wireless is off. Pairing is not allowed.	Auto-Send Not Complete Transfer still pending. Ensure paired device is within range and turned on.
A meter setting was changed while in Flight Mode .	Pairing to a device cannot be performed while in Flight Mode .	The blood glucose result has not been sent to the default paired
The setting change will not take effect until Flight Mode is	Retry pairing when the meter is not in Flight Mode .	device. The transfer is still pending.
turned off.		Place meter and paired device closer together.

A meter or test strip error has occurred.

Your blood glucose may be extremely high or a meter or a test strip error has occurred.

- If your test result matches how you feel, contact your healthcare professional immediately.
- If your test result does not match how you feel, repeat the blood glucose test. See the Unusual Blood Glucose Results section in the chapter Blood Glucose Tests.
 - If the E-3 code still appears for your blood glucose test, your blood glucose result may be extremely high and above the system's reading range. Contact your healthcare professional immediately.

- If the second test result does not match how you feel, perform a control test with the control solution and a new test strip.
 - If the control result is within the acceptable range, review the proper testing procedure and repeat the blood glucose test with a new test strip.
 - If the control result is not within the acceptable range, see the Understanding Out-of-Range Control Results section in the chapter Control Tests.

8

Electronic Error Remove batteries, press any button, reinsert batteries. (F-7)

An electronic error occurred, or in rare cases, a used test strip was removed and reinserted.

Remove the batteries, press any button, and reinsert the batteries. Perform a blood glucose or control test.

Temperature Error Move meter to correct temperature and wait 5 minutes. (E-8)

The temperature is above or below the proper range for the system.

Refer to the test strip package insert for system operating conditions. Move to an area with the appropriate conditions, wait 5 minutes, and repeat the blood glucose or control test. Do not artificially heat or cool the meter.

Change the batteries now. If the message reappears after the batteries have been replaced, slide the battery drawer out of the meter, press any meter button, then reinsert the battery drawer into position.

(x) Time/Date	
Error	
Correct	
time/date	
if necessary.	
(E-10)	
The time and date setting may be inc	correct.
Make sure the time and date are cor adjust, if necessary.	rect and

9

See the literature packaged with the test strips and control solutions for the latest information on product specifications and limitations.

Specifications	
Blood volume	Refer to the test strip package insert.
Sample type	
Measuring time	
Measurement range	
Test strip storage conditions	
System operating conditions	
Meter storage conditions	Temperature: -25–70 °C
Memory capacity	750 blood glucose results and 30 control results with time and date
Automatic off	90 seconds
Power supply	Two 3-volt lithium batteries (coin cell type CR2032)
Display	LCD
Dimensions	$80 \times 47 \times 19$ mm (LWH)
Weight	Approx. 40 g (with batteries)
Construction	Hand-held
Protection class	III
Meter type	The Accu-Chek Aviva Connect meter is suitable for continuous operation.
Control solution storage conditions	Refer to the control solution package insert.

Bluetooth® Wireless Technology

This device complies with Industry Canada license-exempt RSS standard(s). Operation is subject to the following two conditions: (1) this device may not cause interference, and (2) this device must accept any interference, including interference that may cause undesired operation of the device.

Under Industry Canada regulations, this radio transmitter may only operate using an antenna of a type and maximum (or lesser) gain approved for the transmitter by Industry Canada. To reduce potential radio interference to other users, the antenna type and its gain should be so chosen that the equivalent isotropically radiated power (e.i.r.p.) is not more than that necessary for successful communication.

Industry Canada ICES-003 Compliance Label:

CAN ICES-3 (B)/NMB-3(B)

Electromagnetic Compatibility – This meter meets the electromagnetic immunity requirements as per EN ISO 15197. The chosen basis for electrostatic discharge immunity testing was basic standard IEC 61000-4-2. In addition, the meter meets the electromagnetic emissions requirements as per EN 61326. The meter's electromagnetic emission is thus low. Interference from the meter to other electrically-driven equipment is not anticipated.

Performance Analysis - Refer to the test strip package insert.

Test Principle - Refer to the test strip package insert.

Declaration of Conformity – Roche Diagnostics hereby declares that the Accu-Chek Aviva Connect blood glucose meter conforms with the basic requirements and other relevant regulations of the European Directive 1999/5/EC. The conformity declaration may be found at the following Internet address: http://declarations.accu-chek.com

Product Safety Information

\land WARNING

- Strong electromagnetic fields may interfere with the proper operation of the meter. Do not use the meter close to sources of strong electromagnetic radiation.
- To avoid electrostatic discharge, do not use the meter in a very dry environment, especially one in which synthetic materials are present.

9

A WARNING

- During blood glucose testing, the meter itself may come into contact with blood. Used meters therefore
 carry a risk of infection. Before discarding the meter, remove the battery or batteries. Discard used
 meters according to the regulations applicable in your country. Contact the local council and authority
 for information about correct disposal.
- The meter falls outside the scope of the European Directive 2002/96/EC Directive on waste electrical and electronic equipment (WEEE).
- Discard used batteries according to local environmental regulations.

Explanation of Symbols

9

These symbols may appear on the packaging, on the type plate, and in the instructions for the Accu-Chek Aviva Connect meter.

Ĩ	Consult instructions for use
\triangle	Caution, refer to safety-related notes in the instructions for use accompanying this product.
1	Temperature limitation (store at)
	Manufacturer
REF	Catalogue number
IVD	In vitro diagnostic medical device
C E 0088	This product fulfils the requirements of the European Directive 98/79/EC on in vitro diagnostic medical devices.
•	3-volt coin cell type CR2032

Additional Supplies

Test Strips: Accu-Chek Aviva test strips

Control Solutions: Accu-Chek Aviva control solutions

9

Low Blood Glucose Index or High Blood Glucose Index^{1, 2, 3}

These figures represent the frequency and the resulting risk of blood glucose values being too low or too high. Figures should be as low as possible.

The following table provides an overview to assess the risk of blood glucose values being too low or too high:

Risk	Low blood glucose index	High blood glucose index
minimal	≤1,1	≤5,0
low	1,1–2,5	5,0–10,0
medium	2,5–5,0	10,0–15,0
high	>5,0	>15,0

A WARNING

The index values for low blood glucose or high blood glucose in the table are not blood glucose values. Ask your healthcare professional if you want to change your therapy based on the index values.

References

¹ Boris P. Kovatchev, Martin Straume, Daniel J. Cox, Leon S. Farhy (2001) "Risk analysis of blood glucose data: a quantitative approach to optimizing the control of insulin dependent diabetes." *Journal of Theoretical Medicine*, 3: pp 1-10.

² Boris P. Kovatchev, Daniel J. Cox, Anand Kumar, Linda Gonder-Frederick, William L. Clarke (2003) "Algorithmic Evaluation of Metabolic Control and Risk of Severe Hypoglycemia in Type 1 and Type 2 Diabetes Using Self-Monitoring Blood Glucose Data." *Diabetes Technology & Therapeutics*, 5(5): pp 817-828.

^a Boris P. Kovatchev (2006) "Is Glycemic Variability Important to Assessing Antidiabetes Therapies." *Current Diabetes Reports*, 6: pp 350-356.

Guarantee

The statutory provisions on rights in consumer goods sales in the country of purchase shall apply.

Guarantee

Index

A

alternate site testing 26 averages 44

B

batteries, changing 77 battery type 87 blood glucose index 91 blood glucose results, unusual 27 blood glucose test, performing 20 button, functions 7

C

control solution 13 control test, performing 14

D

data transfer default 62

G

guarantee 93

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high blood glucose 28 hyperglycemia 28 hypoglycemia 28

I

icons 10 icons, blood glucose result 24

L

logbook 42 low blood glucose 28

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maintenance, meter 77 meter, cleaning 78

Ρ

pc reports default 63 list report 74 standard day report 72 standard week report 73 trend report 71 product limitations 87 product safety information 88 product specifications 87

S

settings, meter 29 beeper 32 language 39 last result 38 markers 35 target range 35 time/date 31 USB default 34 symbols 90

T

technical information 87 test strips 20

U

use by date 13, 20

V

view data, pc internet browser 65 view data in software 64 view pc reports 67

W

wireless communication 47 auto-send 52 default device 51 delete paired device 57 flight mode 50 pair additional device 54 pairing 48 sync time 53

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ACCU-CHEK, ACCU-CHEK AVIVA, and ACCU-CHEK AVIVA CONNECT are trademarks of Roche. ACCU-CHEK, ACCU-CHEK AVIVA et ACCU-CHEK AVIVA CONNECT sont des marques de Roche.
The *Bluetooth*[®] word mark and logos are registered trademarks owned by Bluetooth SIG, Inc. and any use of such marks by Roche is under license. Other trademarks and trade names are those of their respective owners.
Roche Diagnostics GmbH Sandhofer Strasse 116 68305 Mannheim, Germany www.accu-chek.com